**Infrarot Therapielampe**  
**Infrared therapy lamp**  
**Lampe de thérapie infrarouge**  
**Infrarood therapielamp**  
**Lampa do terapii na podczerwień**  
**Lámpara de terapia de infrarrojos**

**IL-100**

DE

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EN

## Instruction Manual Infrared therapy lamp

Voltage: 230V~ 50Hz  
Power: 100W

## **BENEFITS OF INFRARED HEAT**

The beneficial effect of infrared heat is well known. It actually penetrates below the skin to stimulate blood flow which helps to speed heating, relieve minor muscular pain, relax muscles, and reduce joint stiffness. Recommended in cases such as backache, tennis elbow, muscle strains and soreness. Indispensable for athletes, gardeners, busy Moms, and active people of all ages.

## **USAGE**

1. The infrared lamp should always be on a stable surface, such as a table. The infrared heater should not be aimed at textiles. The body parts to be irradiated should be free of metal objects. Belt buckles, bra buckles and piercings should be removed beforehand.
2. Connect the infrared lamp to the power supply, switch on the lamp and relax in front of the device.
3. The maximum duration of treatment should not exceed 15 minutes per session. Observe the recommended distance of 30-40 cm.
4. For the first application, it is advisable to reduce the treatment time to eight to ten minutes. This makes it easier for the skin to adjust to the infrared radiation.
5. Point the infrared lamp at the painful part of the body (e.g. back, shoulder, face).
6. Close or cover your eyes during treatment.
7. Relax during the application, for example with a quiet CD.

## **IMPORTANT SAFETY MEASURES**

1. To protect against electric shock, do not use while bathing.
2. Do not immerse the device in water.
3. This device must be kept dry and must not be operated when it is wet or damp.
4. Do not operate if the device or the power cord are damaged.
5. Never pull on the cable, never use the cable as a handle and avoid kinking or twisting the cable. Immediately discontinue use if damage is visible.
6. To avoid overheating and the possibility of burns or fire, use the device only as directed in these instructions and precautions and never fall.

asleep while using the device. Never wrap duvets or other materials around the device.

## **ATTENTION**

1. Temperature sufficiently high that burns can occur.
2. Prevent the eyes from being irradiated.
3. Do not use on pregnant women or children.
4. Do not use for more than 15 minutes at a time.
5. Observe the recommended distance of 30-40cm.
6. The device may only be dismantled and repaired by qualified specialists.

## **CLEANING**

Regularly remove dust with a dry cloth.

## Waste Management

The symbol above and on the product means that the product is classed as Electrical or Electronic equipment and should not be disposed with other household or commercial

waste at the end of its useful life.

The Waste of Electrical and Electronic Equipment (WEEE) Directive (2012/19/EU) has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill. Contact local authorities for information on the correct disposal of Electrical or Electronic equipment.

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